

**AdipoGen**<sup>®</sup>   
LIFE SCIENCES

## MANUAL

# Clusterin (human) Competitive ELISA Kit

*For research use only. Not for diagnostic use.*

Version 3 (05-May-2015)

**Cat. No. AG-45A-0013YEK-KI01**

[www.adipogen.com](http://www.adipogen.com)

## Table of Contents

---

<b>1. Intended Use</b>	<b>3</b>
<b>2. Introduction</b>	<b>3</b>
<b>3. General References</b>	<b>4</b>
<b>4. Assay Principle</b>	<b>5</b>
<b>5. Handling &amp; Storage</b>	<b>5</b>
<b>6. Kit Components</b>	<b>5</b>
<b>7. Materials Required but <i>Not</i> Supplied</b>	<b>6</b>
<b>8. General ELISA Protocol</b>	<b>7</b>
8.1. Preparation and Storage of Reagents	7
8.2. Sample Collection, Storage and Dilution	8
8.3. Assay Procedure (Checklist)	9
<b>9. Calculation of Results</b>	<b>10</b>
<b>10. Typical Data</b>	<b>10</b>
<b>11. Performance Characteristics</b>	<b>11-12</b>
<b>12. Technical Hints and Limitations</b>	<b>13</b>
<b>13. Troubleshooting</b>	<b>14</b>
<b>14. Notes</b>	<b>15</b>
<hr/>	
<b>Product Specific References</b>	<b>Backcover</b>

---

## 1. Intended Use

The Clusterin (human) Competitive ELISA Kit is to be used for the *in vitro* quantitative determination of human clusterin in serum, plasma and cell culture supernatant. This ELISA Kit is for research use only.

## 2. Introduction

Clusterin is a multifunctional protein, which has many alternative names; apolipoprotein J (ApoJ), sulfated glycoprotein 2 (SGP2), and complement-associated protein SP-40. Apolipoprotein J is the human analog of the rat protein present in high concentrations in the testis, sulfated glycoprotein-2. It is a 70-kD protein associated with high-density lipoproteins (HDL) in human plasma. Its primary structure was deduced by de Silva et al. (1) using the combined strategies of protein sequencing and cDNA cloning and sequencing. There is a single copy of the APOJ gene in the human and mouse genomes. The protein is synthesized as a 427-amino acid polypeptide that is posttranslationally cleaved at an internal bond between arg205 and ser206. Two subunits, designated alpha (34 to 36 kD), corresponding to residues 1-205, and beta (36 to 39 kD), corresponding to residues 206-427, are associated through disulfide bonds. Studies indicated that the alpha and beta subunits are derived from a common precursor by proteolytic cleavage and that the subunits, while distinct, have limited regions of homology. De Silva et al. (2) found APOJ mRNA (1.9 kb) in all but one tissue examined. Its concentration was relatively high in brain, ovary, testis, and liver, lower in heart, spleen, lung, and breast, and absent in T lymphocytes. Apolipoprotein J is distinct from other known apolipoproteins in molecular weight, subunit structure, and isoelectric point. Clusterin is a major serum protein whose normal concentrations are around 100 ug/ml. Clusterin is a member of the human complement system by directly demonstrating its presence within the S-protein-containing soluble variant of the C5b-9 complex, SC5b-9. It acts as a control mechanism of the complement cascade; specifically, it prevents the binding of a C5b-C7 complex to the membrane of the target cell and in this way inhibits complement-mediated cytotoxicity. The findings of Kirszbaum et al. (3) document a link between the immune and reproductive systems. For this reason the term clusterin is used for the protein in both human serum and seminal fluid. Clusterin is encoded by a single gene, but there are two distinct forms of clusterin; 1) nuclear form and 2) secretory form. The secretory form exists in both serum and semen. The secretory form undergoes a heavy protein processing including cleavage and glycosylation as described above. The nuclear form is a simple open reading frame initiated from the 2nd initiation codon whereas the secretory form initiates from the 1st initiation codon in the leader peptide. While the secretory form of clusterin is cytoprotective, the nuclear form is apoptotic. The serum levels of the secretory form become increasing due to many cancers (4). Clusterin is related to neuron. There are a growing number of papers regarding the relationship between Alzheimer Disease (AD) and clusterin. Clusterin acts on microglial cells, which are brain macrophage-equivalent (4). Thus, clusterin could be a novel diagnostic marker in many major human diseases like cancers or neurodegenerative diseases.

### 3. General References

- (1) Apolipoprotein J: structure and tissue distribution: H.V. de Silva, et al.; *Biochemistry* **29**, 5380 (1990)
- (2) Purification and characterization of apolipoprotein J: H.V. de Silva, et al.; *J. Biol. Chem.* **265**, 14292 (1990)
- (3) Molecular cloning and characterization of the novel, human complement-associated protein, SP-40,40: a link between the complement and reproductive systems: L. Kirschbaum, et al.; *EMBO J.* **8**, 711 (1989)
- (4) Challenge and promise: roles for clusterin in pathogenesis, progression and therapy of cancer: B. Shannan, et al.; *Cell Death & Differentiation* **13**, 12 (2006)

## 4. Assay Principle

This assay is a competitive Enzyme Linked-Immunosorbent Assay (ELISA) for quantitative determination of human clusterin in biological fluids. A polyclonal antibody recognizing native human clusterin reacts with a series of predetermined recombinant human clusterin standard proteins or samples under competition in the human clusterin-coated plate. Their relative reactivity is plotted with that of the standard proteins.

## 5. Handling & Storage

- Reagent must be stored at 2-8°C when not in use.
- Plate and reagents should be at room temperature before use.
- Do not expose reagents to temperatures greater than 25°C.

## 6. Kit Components

1 plate coated with clusterin (h) Recombinant protein	(6 x 16-well strips)	
2 bottles Wash Buffer 10X	(2 x 30 ml)	(Wash Buffer 10X)
2 bottles ELISA Buffer 10X	(2 x 30 ml)	(ELISA Buffer 10X)
1 vial Detection Antibody	(20 µl)	(DET)
1 vial HRP 100X (HRP Conjugated anti-rabbit IgG)	(150 µl)	(HRP 100X)
1 vial human clusterin Standard (lyophilized)	(5 µg)	(STD)
1 bottle TMB Substrate Solution	(12 ml)	(TMB)
1 bottle Stop Solution	(12 ml)	(STOP)
2 plate sealers (plastic film)		
2 silica Gel Minibags		

## 7. Materials Required but *Not* Supplied

- Microtiterplate reader at 450 nm
- Calibrated precision single and multi-channel pipettes. Disposable pipette tips
- Deionized water
- Microtubes or equivalent for preparing dilutions
- Disposable plastic containers for preparing working buffers
- Plate washer: automated or manual
- Glass or plastic tubes for diluting and aliquoting standard

## 8. General ELISA Protocol

### 8.1. Preparation and Storage of Reagents

**NOTE:** Prepare just the appropriate amount of the buffers necessary for the assay.

- **Wash Buffer 10X** has to be diluted with deionized water 1:10 before use (e.g. 50 ml Wash Buffer 10X + 450 ml water) to obtain Wash Buffer 1X.
- **ELISA Buffer 10X** has to be diluted with deionized water 1:10 before use (e.g. 20 ml ELISA Buffer 10X + 180 ml water) to obtain ELISA Buffer 1X.
- **Detection Antibody (DET)** has to be diluted to 1:1660 in ELISA Buffer 1X (6 µl DET + 10 ml ELISA Buffer 1X).

**NOTE:** The diluted Detection Antibody is not stable and cannot be stored!

- **HRP 100X (HRP Conjugated anti-rabbit IgG)** has to be diluted to the working concentration by adding 100 µl in 10 ml of ELISA Buffer 1X (1:100).

**NOTE:** The diluted HRP is used within one hour of preparation.

- **Human clusterin Standard (STD)** has to be reconstituted with 1 ml of deionized water.
  - This reconstitution produces a stock solution of 5 µg/ml. Mix the standard to ensure complete reconstitution and allow the standard to sit for a minimum of 15 minutes. Mix well prior to making dilutions.

**NOTE:** The reconstituted standard is aliquoted and stored at -20°C.

- Dilute the standard protein concentrate (STD) (5 µg/ml) in ELISA Buffer 1X. A seven-point standard curve in ELISA Buffer 1X is recommended.
- Suggested standard points are:  
**5, 2.5, 1, 0.5, 0.25, 0.1, 0.01 and 0.001 µg/ml.**

**Dilute further for the standard curve:**

To obtain	Add	Into
5 µg/ml	-	-
2.5 µg/ml	300 µl of clusterin (5 µg/ml)	300 µl of ELISA Buffer 1X
1 µg/ml	200 µl of clusterin (2.5 µg/ml)	300 µl of ELISA Buffer 1X
0.5 µg/ml	300 µl of clusterin (1 µg/ml)	300 µl of ELISA Buffer 1X
0.25 µg/ml	300 µl of clusterin (0.5 µg/ml)	300 µl of ELISA Buffer 1X
0.1 µg/ml	200 µl of clusterin (0.25 µg/ml)	300 µl of ELISA Buffer 1X
0.01 µg/ml	50 µl of clusterin (0.1 µg/ml)	450 µl of ELISA Buffer 1X
0.001 µg/ml	50 µl of clusterin (0.01 µg/ml)	450 µl of ELISA Buffer 1X

## 8.2. Sample Collection, Storage and Dilution

**Serum** : Use a serum separator tube. Let samples clot at room temperature for 30 minutes before centrifugation for 20 minutes at 1,000xg. Assay freshly prepared serum or store serum in aliquot at ≤ -20°C for later use. Avoid repeated freeze/thaw cycles.

**Plasma** : Collect plasma using heparin, EDTA, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000xg within 30 minutes of collection. Assay freshly prepared plasma or store plasma sample in aliquot at ≤ -20°C for later use. Avoid repeated freeze/ thaw cycles.

**Serum, Plasma or Cell Culture Supernatant** have to be diluted in ELISA Buffer 1X. Samples containing visible precipitates must be clarified before use.

**NOTE:** As a starting point, 1/500 dilution of serum or plasma is recommended! If sample values fall outside the detection range of the assay, a lower or higher dilution may be required!



### 8.3. Assay Procedure (Checklist)


<input type="checkbox"/>	<p>1. Determine the number of 16-well strips needed for the assay and insert them in the frame for current use. The extra strips should be resealed in the foil pouch bag and stored at 4°C.</p> <p><b>NOTE:</b> Remaining 16-well strips coated with clusterin protein when opened can be stored at 4°C for up to 1 month.</p>
<input type="checkbox"/>	<p>2. Add 50 µl of the different standards into the appropriate wells in duplicate! At the same time, add 50 µl of diluted serum, plasma or cell culture supernatant samples in duplicate to the wells (<b>see 8.1. Preparation and Storage of Reagents and 8.2. Preparation of Samples</b>).</p>
<input type="checkbox"/>	<p>3. Add 50 µl to each well of the Detection Antibody (<b>DET</b>) (<b>see 8.1. Preparation and Storage of Reagents</b>) and tap gently on the side of the plate to mix.</p>
<input type="checkbox"/>	<p>4. Cover the plate with plate sealer and incubate for <b>1 hour at 37°C</b>.</p>
<input type="checkbox"/>	<p>5. Aspirate the coated wells and add 300 µl of Wash Buffer 1X using a multi-channel pipette or auto-washer. Repeat the process for a total of three washes. After the last wash, complete removal of liquid is essential for good performance.</p>
<input type="checkbox"/>	<p>6. Add 100 µl to each well of the diluted HRP Conjugated anti-rabbit IgG (<b>HRP</b>) (<b>see 8.1. Preparation and Storage of Reagents</b>).</p>
<input type="checkbox"/>	<p>7. Cover the plate with plate sealer and incubate for <b>1 hour at 37°C</b>.</p>
<input type="checkbox"/>	<p>8. Aspirate the coated wells and add 300 µl of Wash Buffer 1X using a multi-channel pipette or auto-washer. Repeat the process for a total of five washes. After the last wash, complete removal of liquid is essential for good performance.</p>
<input type="checkbox"/>	<p>9. Add 100 µl to each well of TMB Substrate Solution (<b>TMB</b>).</p>
<input type="checkbox"/>	<p>10. Allow the color reaction to develop <b>at room temperature (RT°C) in the dark for 20 minutes</b>.</p>
<input type="checkbox"/>	<p>11. Stop the reaction by adding 100 µl of Stop Solution (<b>STOP</b>). Tap the plate gently to ensure thorough mixing. The substrate reaction yields a blue solution that turns yellow when Stop Solution (<b>STOP</b>) is added.</p>
<input type="checkbox"/>	<p><b>! CAUTION: CORROSIVE SOLUTION!</b></p>
<input type="checkbox"/>	<p>12. Measure the OD at 450 nm in an ELISA reader within 30 minutes.</p>

## 9. Calculation of Results

- Average the duplicate readings for each standard, control and sample.
- Generate the standard curve by plotting the average absorbance obtained for each standard concentration on the vertical (Y) axis vs. the corresponding clusterin concentration (µg/ml) on the horizontal (X) axis (see 10. TYPICAL DATA).
- Calculate the clusterin concentrations of samples by interpolation of the regression curve formula as shown above in a form of a 4-parameter logistic equation.
- If the test samples were diluted, multiply the interpolated values by the dilution factor to calculate the concentration of human clusterin in the samples.

## 10. Typical Data

The following data are obtained using the different concentrations of standard as described in this protocol:



Standard hClusterin (µg/ml)	Optical Density (mean)
5	0.550
2.5	0.631
1	0.896
0.5	1.179
0.25	1.437
0.1	1.851
0.01	2.678
0.001	2.845

**Figure:** Standard curve

## 11. Performance Characteristics

### A. Sensitivity (Limit of detection):

The lowest level of clusterin that can be detected by this assay is 1 ng/ml. **NOTE:** *The Limit of detection was measured by adding two standard deviations to the mean value of 50 zero standard.*

**B. Assay range:** 0.001 µg/ml – 5 µg/ml

### C. Specificity:

This ELISA is specific for the measurement of natural and recombinant human clusterin. It does not cross-react with mouse clusterin, human adiponectin, human resistin, human RBP4, human RELM-β, human FABP4, human GPX3, human Nampt, human PAI-1, human leptin, human IL-23, human TNF-α, human ANGPTL3, human ANGPTL4, human ANGPTL7.

### D. Intra-assay precision:

Five samples of known concentrations of human clusterin were assayed in replicates 8 times to test precision within an assay.

Samples	Means (µg/ml)	SD	CV (%)	n
1	130.61	12.83	9.82	8
2	61.72	3.13	5.07	8
3	66.87	4.68	7.02	8
4	70.35	5.31	7.55	8
5	35.94	3.26	9.07	8

### E. Inter-assay precision:

Five samples of known concentrations of human clusterin were assayed in 8 separate assays to test precision between assays.

Samples	Means (µg/ml)	SD	CV (%)	n
1	68.33	6.67	9.76	8
2	51.82	3.47	6.70	8
3	83.02	8.17	9.84	8
4	63.92	5.88	9.20	8
5	49.22	3.53	7.18	8

**F. Linearity:**

Different human serum samples containing clusterin were diluted several fold (1/500 to 1/2,000) and the measured recoveries ranged from 98% to 120%.

Samples	Sample Dilution	Expected (µg/ml)	Observed (µg/ml)	% of Expected
1	1 : 500	153.77	153.77	100
	1 : 1,000	76.89	79.09	103
	1 : 2,000	38.44	37.61	98
2	1 : 500	165.26	165.26	100
	1 : 1,000	82.63	97.13	118
	1 : 2,000	41.32	48.63	118
3	1 : 500	83.24	83.24	100
	1 : 1,000	41.62	45.02	108
	1 : 2,000	20.81	22.48	108

**G. Expected values:**

Clusterin levels range in plasma and serum from **30 to > 200 µg/ml** (from healthy donors).


## 12. Technical Hints and Limitations

- It is recommended that all standards, controls and samples be run in duplicate.
- Do not combine leftover reagents with those reserved for additional wells.
- Reagents from the kit with a volume less than 100 µl should be centrifuged.
- Residual wash liquid should be drained from the wells after last wash by tapping the plate on absorbent paper.
- Crystals could appear in the 10X solution due to high salt concentration in the stock solutions. Crystals are readily dissolved at room temperature or at 37°C before dilution of the buffer solutions.
- Once reagents have been added to the 16-well strips, DO NOT let the strips DRY at any time during the assay.
- Keep TMB Substrate Solution (TMB) protected from light.
- The Stop Solution (STOP) consists of sulfuric acid. Although diluted, the Stop Solution (STOP) should be handled with gloves, eye protection and protective clothing.

## 13. Troubleshooting

PROBLEM	POSSIBLE CAUSES	SOLUTIONS
No signal or weak signal	Omission of key reagent	Check that all reagents have been added in the correct order.
	Washes too stringent	Use an automated plate washer if possible.
	Incubation times inadequate	Incubation times should be followed as indicated in the manual.
	Plate reader settings not optimal	Verify the wavelength and filter setting in the plate reader.
	Incorrect assay temperature	Use recommended incubation temperature. Bring substrates to room temperature before use.
High background	Concentration of HRP too high	Use recommended dilution factor.
	Inadequate washing	Ensure all wells are filling wash buffer and are aspirated completely.
Poor standard curve	Wells not completely aspirated	Completely aspirate wells between steps.
	Reagents poorly mixed	Be sure that reagents are thoroughly mixed.
Unexpected results	Omission of reagents	Be sure that reagents were prepared correctly and added in the correct order.
	Dilution error	Check pipetting technique and double-check calculations.

## 14. Notes



#### **Product Specific References:**

1. R.T. Blankley, et al.; BJOG. **116**, 1473 (2009)
2. H.R. Yu, et al.; Pediatr. Allergy Immunol. **20**, 699 (2009)
3. S. Pucci, et al.; Am. J. Gastroenterol. **104**, 2807 (2009)
4. H.R. Yu, et al.; Pediatr. Cardiol. **31**, 1151 (2010)
5. D. Kevans, et al.; J. Gastrointest. Cancer **43**, 305 (2012)
6. J. Magalhaes & M.J. Saraiva; Amyloid **19**, 3 (2012)
7. J.C. Won, et al.; PLOS One **9**, e103351 (2014)

*For more References please visit [www.adipogen.com](http://www.adipogen.com)!*

**Adipogen Life Sciences**  
Schützenstrasse 12  
CH-1410 Liestal  
Switzerland  
TEL: +41-61-926-60-40  
FAX: +41-61-926-60-49  
Email: [info@adipogen.com](mailto:info@adipogen.com)

**AdipoGen**<sup>®</sup>   

---

**LIFE SCIENCES**