

**AdipoGen**<sup>TM</sup> 

## **MANUAL**

### **CTRP5 (human) Competitive ELISA Kit**

*For research use only. Not for diagnostic use.*

Version 2 (14-March-2011)

**Cat. No. AG-45A-0031EK-KI01**

[www.adipogen.com](http://www.adipogen.com)

## Table of Contents

---

<b>1. Intended Use</b>	<b>3</b>
<b>2. Introduction</b>	<b>3</b>
<b>3. General References</b>	<b>4</b>
<b>4. Assay Principle</b>	<b>5</b>
<b>5. Handling &amp; Storage</b>	<b>5</b>
<b>6. Kit Components</b>	<b>5</b>
<b>7. Materials Required but <i>Not</i> Supplied</b>	<b>6</b>
<b>8. General ELISA Protocol</b>	<b>7</b>
8.1. Preparation and Storage of Reagents	7
8.2. Sample Collection, Storage and Dilution	8
8.3. Assay Procedure (Checklist)	9
<b>9. Calculation of Results</b>	<b>10</b>
<b>10. Typical Data</b>	<b>10</b>
<b>11. Performance Characteristics</b>	<b>11-12</b>
<b>12. Technical Hints and Limitations</b>	<b>13</b>
<b>13. Troubleshooting</b>	<b>14</b>
<b>14. Assay Flow Chart</b>	<b>15</b>
<b>Product Specific References</b>	<b>Backcover</b>

---

## 1. Intended Use

The CTRP5 (human) Competitive ELISA Kit is to be used for the *in vitro* quantitative determination of human CTRP5 in serum, plasma and cell culture supernatant. This ELISA Kit is for research use only.

## 2. Introduction

CTRP5 is a member of the C1q and tumor necrosis factor superfamily whose structure resembles adiponectin, being a 25 kD secretory protein (1-2). CTRP5 was identified in a cDNA library enriched for genes showing expression specific for the retinal pigment epithelium (RPE) (3). A subsequent RT-PCR study demonstrated CTRP5 expression in RPE, liver, lung, placenta, and brain. It has been proposed that a CTRP mutation (S163R) plays a critical role in affecting its higher order protein structure, potentially leading to a cause of abnormal adhesion between the RPE and Bruch membrane (4-5). A recent data showed that CTRP is one of the genes highly induced by elimination of mitochondria and was able to activate AMPK in a rat myotube cell line, L6. Stimulation of L6 with recombinant CTRP5, full length as well as globular domain, enhanced glucose uptake and fatty acid oxidation. These biochemical events did not seem to be mediated via AdipoR1 or AdipoR2, suggesting that a novel receptor(s) may exist for CTRP5 in this muscle cell line (6). Some CTRP members can physically interact with adiponectin, forming various multimeric structures (2). A compensatory upregulation of mouse CTRP members was observed in adiponectin-deficient mice (1). Therefore, measuring serum or plasma CTRP5 may provide important information on its involvement in novel metabolism.

### 3. General References

- (1) Molecular, biochemical and functional characterizations of C1q/TNF family members: adipose-tissue-selective expression patterns, regulation by PPAR-gamma agonist, cysteine-mediated oligomerizations, combinatorial associations and metabolic functions: G.W. Wong, et al.; *Biochem. J.* 416, 161 (2008)
- (2) Identification and characterization of CTRP9, a novel secreted glycoprotein, from adipose tissue that reduces serum glucose in mice and forms heterotrimers with adiponectin: G.W. Wong, et al.; *FASEB J.* 23, 241 (2009)
- (3) CTRP5 is a membrane-associated and secretory protein in the RPE and ciliary body and the S163R mutation of CTRP5 impairs its secretion: M.N. Mandal, et al.; *Invest. Ophthalmol. Vis. Sci.* 47, 5505 (2006)
- (4) Spatial and temporal expression of MFRP and its interaction with CTRP5: M.N. Mandal, et al.; *Invest. Ophthalmol. Vis. Sci.* 47, 5514 (2006)
- (5) Mutation in a short-chain collagen gene, CTRP5, results in extracellular deposit formation in late-onset retinal degeneration: a genetic model for age-related macular degeneration: C. Hayward, et al.; *Hum. Mol. Genet.* 12, 2657 (2003)
- (6) C1q tumor necrosis factor alpha-related protein isoform 5 is increased in mitochondrial DNA-depleted myocytes and activates AMP-activated protein kinase: S.Y. Park, et al.; *J. Biol. Chem.* 284, 27780 (2009)

## 4. Assay Principle

This assay is a competitive Enzyme Linked-Immunosorbent Assay (ELISA) for quantitative determination of human CTRP5 in biological fluids. A polyclonal antibody recognizing native human CTRP5 reacts with a series of predetermined recombinant human CTRP5 standard proteins or samples under competition in the human CTRP5-coated plate. Their relative reactivity is plotted with that of the standard proteins.

## 5. Handling & Storage

- Reagent must be stored at 2-8°C when not in use.
- Plate and reagents should be at room temperature before use.
- Do not expose reagents to temperatures greater than 25°C.

## 6. Kit Components

1 plate coated with human CTRP5 Recombinant Protein	(12 x 8-well strips)
1 bottle Wash Buffer 10X	(50 ml)
1 bottle Diluent 5X	(50 ml)
1 bottle Detection Antibody	(12 ml)
1 vial Detector 100X (HRP Conjugated anti-rabbit IgG)	(150 µl)
1 vial human CTRP5 Standard (lyophilized)	(5 µg)
1 vial human CTRP5 QC sample (lyophilized)	
1 bottle TMB Substrate Solution	(12 ml)
1 bottle Stop Solution	(12 ml)
3 plate sealers (plastic film)	

## 7. Materials Required but *Not* Supplied

- Microtiterplate reader at 450 nm, with the correction wavelength set at 540 nm or 570 nm
- Calibrated precision single and multi-channel pipettes. Disposable pipette tips
- Deionized water
- Microtubes or equivalent for preparing dilutions
- Disposable plastic containers for preparing working buffers
- Plate washer: automated or manual
- Glass or plastic tubes for diluting and aliquoting standard

## 8. General ELISA Protocol

### 8.1. Preparation and Storage of Reagents

**NOTE:** Prepare just the appropriate amount of the buffers necessary for the assay.

- **Wash Buffer 10X** has to be diluted with deionized water 1:10 before use (e.g. 50 ml Wash Buffer 10X + 450 ml water) to obtain Wash Buffer 1X.
- **Diluent 5X** has to be diluted with deionized water 1:5 before use (e.g. 50 ml Diluent 5X + 200 ml water) to obtain Diluent 1X.
- **Detector 100X (HRP Conjugated anti-rabbit IgG)** has to be diluted to the working concentration by adding 120 µl in 12 ml of Diluent 1X (1:100).

**NOTE:** The diluted Detector is used within one hour of preparation.


- **Human CTRP5 Standard (STD)** has to be reconstituted with 1 ml of deionized water.
  - This reconstitution produces a stock solution of 5 µg/ml. Mix the standard to ensure complete reconstitution and allow the standard to sit for a minimum of 15 minutes. Mix well prior to making dilutions.

**NOTE:** The reconstituted standard is aliquoted and stored at -20°C.

- Dilute the standard protein concentrate (STD) (**5 µg/ml**) in Diluent 1X. A seven-point standard curve in Diluent 1X is recommended.
  - Suggested standard points are:  
**5 , 2.5 , 1 , 0.5 , 0.25 , 0.1 , 0.01 and 0.001 µg/ml.**
- **Human CTRP5 QC sample** has to be reconstituted with 1 ml of deionized water.
    - Refer to the Certificate of Analysis for current QC sample concentration. Mix the QC sample to ensure complete reconstitution and allow the QC sample to sit for a minimum of 15 minutes. The reconstituted QC sample is ready to use, do not dilute it.

**Dilute further for the standard curve:**

To obtain	Add	Into
5 µg/ml	-	-
2.5 µg/ml	300 µl of CTRP5 (5 µg/ml)	300 µl of Diluent 1X
1 µg/ml	200 µl of CTRP5 (2.5 µg/ml)	300 µl of Diluent 1X
0.5 µg/ml	300 µl of CTRP5 (1 µg/ml)	300 µl of Diluent 1X
0.25 µg/ml	300 µl of CTRP5 (0.5 µg/ml)	300 µl of Diluent 1X
0.1 µg/ml	200 µl of CTRP5 (0.25 µg/ml)	300 µl of Diluent 1X
0.01 µg/ml	50 µl of CTRP5 (0.1 µg/ml)	450 µl of Diluent 1X
0.001 µg/ml	50 µl of CTRP5 (0.01 µg/ml)	450 µl of Diluent 1X



## 8.2. Sample Collection, Storage and Dilution

**Serum** : Use a serum separator tube. Let samples clot at room temperature for 30 minutes before centrifugation for 20 minutes at 1,000xg. Assay freshly prepared serum or store serum in aliquot at ≤ -20°C for later use. Avoid repeated freeze/thaw cycles.

**Plasma** : Collect plasma using heparin, EDTA, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000xg within 30 minutes of collection. Assay freshly prepared plasma or store plasma sample in aliquot at ≤ -20°C for later use. Avoid repeated freeze/ thaw cycles.

**Serum, Plasma or Cell Culture Supernatant** have to be diluted in Diluent 1X. Samples containing visible precipitates must be clarified before use.

**NOTE:** As a starting point, 1/2 dilution of serum or plasma is recommended! If samples fall the outside range of assay, a lower or higher dilution may be required!



### 8.3. Assay Procedure (Checklist)


<input type="checkbox"/>	<p>1. Determine the number of 8-well strips needed for the assay and insert them in the frame for current use. The extra strips should be resealed in the foil pouch bag and stored at 4°C.</p> <p><b>NOTE:</b> Remaining 8-well strips coated with CTRP5 protein when opened can be stored at 4°C for up to 1 month.</p>
<input type="checkbox"/>	<p>2. Add 50 µl of the different standards and reconstituted QC sample into the appropriate wells in duplicate! At the same time, add 50 µl of diluted serum, plasma or cell culture supernatant samples in duplicate to the wells (<b>see 8.1. Preparation and Storage of Reagents and 8.2. Preparation of Samples</b>).</p>
<input type="checkbox"/>	<p>3. Add 50 µl to each well of the Detection Antibody and tap gently on the side of the plate to mix.</p>
<input type="checkbox"/>	<p>4. Cover the plate with plate sealer and incubate for <b>1 hour at 37°C</b>.</p>
<input type="checkbox"/>	<p>5. Aspirate the coated wells and add 300 µl of Wash Buffer 1X using a multi-channel pipette or auto-washer. Repeat the process for a total of three washes. After the last wash, complete removal of liquid is essential for good performance.</p>
<input type="checkbox"/>	<p>6. Add 100 µl to each well of the diluted Detector (<b>see 8.1. Preparation and Storage of Reagents</b>).</p>
<input type="checkbox"/>	<p>7. Cover the plate with plate sealer and incubate for <b>1 hour at 37°C</b>.</p>
<input type="checkbox"/>	<p>8. Aspirate the coated wells and add 300 µl of Wash Buffer 1X using a multi-channel pipette or auto-washer. Repeat the process for a total of five washes. After the last wash, complete removal of liquid is essential for good performance.</p>
<input type="checkbox"/>	<p>9. Add 100 µl to each well of TMB Substrate Solution.</p>
<input type="checkbox"/>	<p>10. Allow the color reaction to develop <b>at room temperature (RT°C) in the dark for 10 minutes</b>.</p>
<input type="checkbox"/>	<p>11. Stop the reaction by adding 100 µl of Stop Solution. Tap the plate gently to ensure thorough mixing. The substrate reaction yields a blue solution that turns yellow when Stop Solution is added.</p>
	<p><b>! CAUTION: CORROSIVE SOLUTION!</b></p>
<input type="checkbox"/>	<p>12. Measure the OD at 450 nm in an ELISA reader within 30 minutes.</p>

## 9. Calculation of Results

- Average the duplicate readings for each standard, QC and sample.
- Generate the standard curve by plotting the average absorbance obtained for each standard concentration on the vertical (Y) axis vs. the corresponding CTRP5 concentration ( $\mu\text{g/ml}$ ) on the horizontal (X) axis (see **10. TYPICAL DATA**).
- Calculate the CTRP5 concentrations of samples by interpolation of the regression curve formula as shown above in a form of a 4-parameter logistic equation.
- If the test samples were diluted, multiply the interpolated values by the dilution factor to calculate the concentration of human CTRP5 in the samples.

## 10. Typical Data

The following data are obtained using the different concentrations of standard as described in this protocol:



Standard hCTR5 ( $\mu\text{g/ml}$ )	Optical Density (mean)
5	0.277
2.5	0.307
1	0.379
0.5	0.475
0.25	0.672
0.1	0.965
0.01	1.707
0.001	1.935

**Figure:** Standard curve

## 11. Performance Characteristics

### A. Sensitivity (Limit of detection):

The lowest level of CTRP5 that can be detected by this assay is 1 ng/ml. **NOTE:** *The Limit of detection was measured by adding two standard deviations to the mean value of 50 zero standard.*

**B. Assay range:** 0.001 µg/ml – 5 µg/ml

### C. Specificity:

This ELISA is specific for the measurement of natural and recombinant human CTRP5. It does not cross-react with mouse CTRP5, human CTRP2 (globular), human CTRP6, human CTRP7 (globular), human CTRP9 (globular), human CTRP10 (globular), mouse CTRP2 (globular), mouse CTRP9 (globular), human adiponectin, human adiponectin (globular), mouse adiponectin, mouse adiponectin (globular), rat adiponectin, rat adiponectin (globular), human RBP4, human Nampt, human vaspin, human GPX3, human ANGPTL3, human progranulin, human leptin.

Human CTRP5 (globular) and CTRP6 (globular) show weakly 20% cross-reactivity in this assay.

### D. Intra-assay precision:

Four samples of known concentrations of human CTRP5 were assayed in replicates 9 times to test precision within an assay.

Samples	Means (µg/ml)	SD	CV (%)	n
1	0.351	0.02	5.87	9
2	0.293	0.03	10.00	9
3	0.250	0.01	4.05	9
4	0.165	0.01	6.06	9

### E. Inter-assay precision:

Four samples of known concentrations of human CTRP5 were assayed in 3 separate assays to test precision between assays.

Samples	Means (µg/ml)	SD	CV (%)	n
1	0.348	0.02	6.36	3
2	0.279	0.03	9.08	3
3	0.253	0.02	8.11	3
4	0.162	0.01	6.77	3

**F. Linearity:**

Different human serum samples containing CTRP5 were diluted several fold (1 to 1/4) and the measured recoveries ranged from 79% to 130%.

Samples	Sample Dilution	Expected (µg/ml)	Observed (µg/ml)	% of Expected
1	1	0.260	0.260	100
	1 : 2	0.130	0.125	95
	1 : 4	0.065	0.071	108
2	1	0.174	0.174	100
	1 : 2	0.087	0.082	94
	1 : 4	0.044	0.056	127
3	1	0.321	0.321	100
	1 : 2	0.161	0.127	79
	1 : 4	0.080	0.071	88

**G. Comparison of serum samples with plasma samples:**

Different human serum samples containing CTRP5 were compared with plasma samples.

Samples	Serum (µg/ml)	Plasma (µg/ml)		
		Citrate	EDTA	Heparin
1	0.127	0.374	0.343	0.126
2	0.098	0.196	0.346	0.098
3	0.125	0.21	0.504	0.134
4	0.145	0.31	0.41	0.152

**H. Expected values:**

CTRP5 levels range in plasma and serum from **0.05 to > 0.5 µg/ml** (from healthy donors).


## 12. Technical Hints and Limitations


- It is recommended that all standards, QC samples and samples be run in duplicate.
- Do not combine leftover reagents with those reserved for additional wells.
- Reagents from the kit with a volume less than 100 µl should be centrifuged.
- Residual wash liquid should be drained from the wells after last wash by tapping the plate on absorbent paper.
- Crystals could appear in the 10X solution due to high salt concentration in the stock solutions. Crystals are readily dissolved at room temperature or at 37°C before dilution of the buffer solutions.
- Once reagents have been added to the 8-well strips, DO NOT let the strips DRY at any time during the assay.
- Keep TMB Substrate Solution protected from light.
- The Stop Solution consists of phosphoric acid. Although diluted, the Stop Solution should be handled with gloves, eye protection and protective clothing.

### 13. Troubleshooting

PROBLEM	POSSIBLE CAUSES	SOLUTIONS
No signal or weak signal	Omission of key reagent	Check that all reagents have been added in the correct order.
	Washes too stringent	Use an automated plate washer if possible.
	Incubation times inadequate	Incubation times should be followed as indicated in the manual.
	Plate reader settings not optimal	Verify the wavelength and filter setting in the plate reader.
	Incorrect assay temperature	Use recommended incubation temperature. Bring substrates to room temperature before use.
High background	Concentration of detector too high	Use recommended dilution factor.
	Inadequate washing	Ensure all wells are filling wash buffer and are aspirated completely.
Poor standard curve	Wells not completely aspirated	Completely aspirate wells between steps.
	Reagents poorly mixed	Be sure that reagents are thoroughly mixed.
Unexpected results	Omission of reagents	Be sure that reagents were prepared correctly and added in the correct order.
	Dilution error	Check pipetting technique and double-check calculations.

## 14. Assay Flow Chart





**Product Specific References:**

*For more References please visit [www.adipogen.com](http://www.adipogen.com)!*

**Adipogen AG**

Schützenstrasse 12  
CH-1410 Liestal  
Switzerland  
TEL: +41-61-926-60-40  
FAX: +41-61-926-60-49  
Email: [info@adipogen.com](mailto:info@adipogen.com)

**Adipogen, Inc.**

Room 401, Venture Building B, Songdo TechnoPark,  
7-50 Songdo-dong, Yeonsu-gu, Incheon,  
Korea 406-840  
TEL: +82-32-858-1470  
FAX: +82-32-831-1470  
Email: [info-kr@adipogen.com](mailto:info-kr@adipogen.com)